



PROXABRUSH® GO-BETWEENS® CLEANERS

Easy to use on-the-go interdental brush. Improves patient compliance and helps remove biofilm from interproximal surfaces.



• Triangular bristles remove up to 25% more plaque¹



• Bendable Neck

- Nylon coated wire to prevent galvanic shock
- Antibacterial bristles²
- Flexible handle and convenient cap

- Patients non-compliant to flossing
- Removal of food particles and plaque
- Patients appreciating on-the-go prevention and management of gingivitis
- Patients concerned with sensitivity caused by galvanic shock
- Periodontal disease gum health maintenance
- Patients with orthodontics, bridge work, implants, or those with fixed retainers

PROXABRUSH® GO-BETWEENS® CLEANERS BULK BOX

IDBBBK6
Contents:
50 Micro Tight
50 Ultra Tight
50 Tight
50 Moderate
50 Wide
25 Extra Wide
Total of 275 units





- Offers a range of sizes from 0.6 to 2.0mm to meet each patient's unique needs
- Convenient, compact format allows quick and easy dispensing
- Reusable storage container helps to keep the operatory organized
- Refill with bulk bags of 50 and save

- Patients non-compliant to flossing
- Removal of food particles and plaque
- Patients appreciating on-the-go prevention and management of gingivitis
- Patients concerned with sensitivity caused by galvanic shock
- Periodontal disease gum health maintenance
- Patients with orthodontics, bridge work, implants, or those with fixed retainers

PROXABRUSH® GO-BETWEENS® ANGLE CLEANERS

Angled design to effectively reach back teeth. Comes with a convenient travel case for on-the-go cleaning.



• Triangular bristles remove up to 25% more plaque¹




• Angled to effectively reach back teeth

- Nylon coated wire to prevent galvanic shock
- Antibacterial bristles²
- Convenient travel case

- Patients who require on-the-go help in the prevention and management of gingivitis
- Patients non-compliant to flossing
- Periodontal disease gum health maintenance and gum massaging on-the-go solution
- Patients concerned with sensitivity caused by galvanic shock
- Patients with orthodontics, bridge work, implants or those with fixed retainers

PROXABRUSH® ORIGINAL SYSTEM

At home original, duel ended handle with secure locking mechanism to hold refills securely.



• Secure lock-in system



• Triangular bristles remove up to 25% more plaque¹

- 2 ended handle for adapting different size refills
- Customized for different patient interdental spaces
- Antibacterial bristles²
- Nylon coated wire to prevent galvanic shock

- Patients who need extended reach to their back teeth
- Patients who have dexterity issues or cannot floss
- Patients requiring a lingual and buccal angled insertion
- Patients who require help in the prevention and management of gingivitis
- Patients concerned with sensitivity caused by galvanic shock
- Patients with orthodontics, bridge work, implants or those with fixed retainers
- Periodontal disease gum health maintenance and gum massaging on-the-go solution
- Patients who need to use different size refills to customize care

SULCUS BRUSH

Dual ended handle and brush helps to remove plaque along the gumline and gingival sulcus.



- Soft antibacterial bristles²
- Two angles allow for easy access to all areas of the mouth including lingual (A) and buccal (B) gumlines
- Easy to use

- Patients who have excess biofilm along gum line
- Patients who have crowns, bridges, braces, lingual wires and implants
- Patients with challenged gumline anatomy due to gum disease, recession or other deficiencies

SOFT-PICKS®

Disposable on-the-go removal of food particles and biofilm.



- Soft rubber bristles safely massage gums
- Tapered design effectively dislodges food and biofilm
- Not made with natural rubber latex

- Patients who require help in the prevention and management of gingivitis
- Patients non-compliant to flossing
- Patients who require on-the go removal of food particles and biofilm
- For patients with orthodontics or fixed retainers, bridge work or implants

1. Compared to conventional round filament. Data on file.
2. Bristles incorporate a patented antibacterial agent for continuous bristle protection during the recommended life of the brush. Bacterial growth that may affect the bristles is inhibited. The agent in the bristles does not protect against disease. As always, rinse your brush after use.